

## Supported Operating Systems:

- Windows 2000 Service Pack 4
- Windows XP Service Pack 3 32 bit
- Windows Vista SP1 32 bit
- Windows 7 Professional SP1 32 Bit (64 bit not tested)
- Windows 7 Ultimate SP1 32 bit and 64 bit
- Windows 7 Ultimate 64-Bit (tested BootCamp 3.0.4 (325) on iMAC 24", Intel Core 2 Duo, 4GB, MAC OS 10.6.7)

Processor	32-bit (x86), min. 1GHz
Main memory	Min. 1GB
Graphic settings	1024x768 pixels (recommended), 16.7 million colors ("TrueColor")
Graphics card	Min. 128MB
Hard disk	At least 5GB free hard disk space for program and 2D data. When using GALILEOS, please take note of the corresponding additional hard disk space required.
USB port	For Version 1.1 and higher, required for USB components only
Browser	Internet Explorer, Version 6.0 SP1 or higher
PDF software	Acrobat Reader, Version >= 8.0 (included on CD). Required for the PDF test report function.

Important: During installation the database has to be connected to the network!

Important: XIO (Intraoral ISA PC plug-in card), XOP (Orthophos ISA PC plug-in card) and IO2 USB (USB module for SIDEXIS intraoral sensors V1.2 and V2.0) are not supported by SIDEXIS XG 2.5.x under 64-bit systems!